



**OLIDENT**  
INNOVATION

# Composite OliRevo

RIGHT SOLUTION FOR EVERY SITUATION



New



## OliREVO


First on the market nano-hybrid composite, that is so convenient and easy to form. Its consistency is called „Wax-like“ enhancing the creation of highly aesthetic restorations with mechanical properties that guarantee the long-term reliable results.

## OliREVolution!


Absolute composite

### LONG-TERM AESTHETIC FILLINGS

Thanks to the exceptional resistance to compression **OliREVO** tolerates high loads, both on the occlusal and incisal surfaces. In addition, **OliREVO** has an unusually low solubility in water - more than seven times lower than the ISO standard! This and the low wear of **OliREVO** provides for long-term esthetic restorations.



OliREVO water solubility index is more than 7 times lower than the ISO standard\*  
\*Manufacturer internal measurements, Germany.



Abrasive wear of OliREVO fillings and competitive materials (200 000 cycles of abrasion ~ approximately one year)\*  
\*Manufacturer internal measurements, Germany.

### ABRASIVE WEAR RESISTANCE

The studies on the fillings surface abrasive wear show that **OliREVO** achieves one of the best results among its competitors. After 200 000 cycles the fillings surface presented 38 $\mu\text{m}$  wear. This proves **OliREVO** restorations are extremely durable and also safe for the opposing teeth.

### „WAX-like“ CONSISTENCY FOR DENTIST COMFORT

The German specialists who created **OliREVO** took great care to make a user friendly composite that will surely impress all users. The unique WAX-like consistency is synonymous with smooth and easy handling of **OliREVO**. The non-stick characteristics allow for efficient molding and carving even in difficult restorative surfaces.

# WAX like

OliREVO - clinical picture by courtesy of Dr. Dominik Morawski, *Denticulus, Nowy Sacz/Limanowa*



Patient with diastema - initial situation.



Diastema closure with OliREVO, tooth 12,11 and 21, 22..